

## Advanced Higher Physics: Assignment Support

### Astronomy & Physics Education Group

### School of Physics & Astronomy

### University of Glasgow

## Measurement of refractive index of water by prism spectrometer

### Introduction – A prism spectrometer

A prism spectrometer is an instrument for observing spectra and measuring angles of deviation of light by a prism. Figure 1 details the key components: collimator, prism table, telescope.




Figure 1: A prism spectrometer

**Collimator:** This is a tube with an adjustable slit at one end and an achromatic converging lens system at the other. The slit should be vertical, and it is usually placed at the focus of the lens so that, when it is illuminated, a beam of parallel light emerges from the collimator. The collimator is fixed to the base of the instrument.

Prism table: This should be horizontal and can be rotated about a vertical axis. A Vernier scale allows the rotation to be measured with respect to the collimator.

Telescope: This is mounted so that it is free to rotate about the same axis as the prism table. It may be focussed to receive parallel light from the collimator. The rotation of the telescope can be measured by another Vernier scale. These are cross-wires in the eyepiece of the telescope, and these should be vertical and horizontal.

The positions of the arms of the spectrometer are measured using angular Vernier scales. A Vernier scale is used when we need to make a measurement of a distance or angle to a greater accuracy than that obtainable though direct visual reading of a linear scale. In this experiment, the Vernier scales mounted on the spectrometer looking something like Figure 2.




Figure 2: A Vernier scale

The position as read from this example is  $38^{\circ} 22'$ . This means “38 degrees 22 arcminutes”. (There are 60 arcminutes in a degree.). How do we get this result?

- The first thing to look for is the position of the ZERO on the Vernier scale. In Figure 3, this is after  $38^{\circ}$ , but before  $39^{\circ}$ .
- To determine position further, look for a point of alignment between the two scales – i.e. a point where the tick marks from both scales light up. In this case, the position is marked by the arrow – at  $22'$  on the Vernier scale.
- The final position is then the addition of these two numbers -  $38^{\circ} 22'$ .

## Measuring $n$ via refraction through a prism

Figure 3 shows one possible arrangement of the prism spectrometer, with light entering the collimator and then being viewed through the telescope. Here the prism has three faces: AB, BC and CA. BC – positioned against the support bracket on the prism table – is ground. Light is refracted at AC and AB. The angle between the incident and emergent direction of the light through the prism is called the angle of deviation,  $\delta$ . If the prism table is rotated clockwise and anti-clockwise the observed spectrum will




Figure 3: Light refracting through a prism

also move. There will come a point, though, where as the prism table is turned the observed spectrum will stop and then move in the opposite direction. The point where this happens is the angle of minimum deviation. When this occurs the rays passing through the prism are parallel to BC.

Once the value of  $\delta$  has been determined, it is possible to calculate the refractive index of the prism material provided the size of the apex angle, A, is known. For the prism used here – an equilateral triangle – it should be  $60^\circ$  but this should be confirmed.

Consider Figure 4, which shows the passage of a monochromatic ray through the prism at minimum deviation. Since the external angle of a triangle equals the sum of the interior opposite angles it follows that

$$\delta = (\theta_a - \theta_g) + (\theta_a - \theta_g) = 2(\theta_a - \theta_g)$$

and

$$A = \theta_g + \theta_g = 2\theta_g \Rightarrow \theta_g = \frac{A}{2}$$

$$\Rightarrow \theta_a = \frac{\delta}{2} + \theta_g = \frac{\delta}{2} + \frac{A}{2} = \frac{\delta + A}{2}$$

Applying Snell's law at the first boundary, and remembering that  $n_a = 1$  we see ...

$$n_g = \frac{n_a \sin(\theta_a)}{\sin(\theta_g)} = \frac{\sin(\theta_a)}{\sin(\theta_g)}$$

$$\Rightarrow n_g = \frac{\sin\left(\frac{\delta + A}{2}\right)}{\sin\left(\frac{A}{2}\right)}$$

[1]




Figure 4: Light passing through prism

# Notes on equipment

## Equipment list

- Spectrometer
- Hollow prism and water source
- Sodium source

## Equipment guidance

### Confirming the apex angle, A

- Turn the prism so that angle A is pointing along the axis of the collimator as shown in Figure 7; make sure that the light is reflected from each of the two faces adjacent to A. (Check it can be seen in the telescope on both sides.) Once happy, fix the prism table in place.
- Note the vernier reading of the telescope on either side of the prism; the difference in these two values will equal  $2A$ .




Figure 7: Confirming A

### Adjusting the position of the prism

- The prism can be rotated by physically turning the prism table by hand – there is also a fine control for making smaller adjustments. Best practice is to use the table to get the prism to roughly where it should be, then switch to the fine control.

### Finding the spectrum

- It can take some time to find the emission spectrum using a prism spectrometer – patience is key. If you are struggling, it is best to aim for an angle of incidence of around  $50^\circ$  for the refraction experiment and around  $45^\circ$  for the reflection

experiment. These don't need to be measured precisely – just gauged by eye to help get started.

## **Focussing the Collimator and Telescope**

### I. Focussing the telescope and cross-wires

- Rotate the telescope to view a point in the lab as far as way as you can.
- Focus the telescope until you get as sharp an image as possible.
- To focus the cross-wires adjust the eyepiece of the telescope until the cross-wires appear sharp. Make sure that the cross-wires are horizontal and vertical. ONCE YOU ARE SURE YOU HAVE THESE FOCUSED, DO NOT TOUCH THE EYEPIECE AGAIN!

### II. Focussing the collimator and telescope using Schuster's method

- Set the equipment so that you can see the line spectrum from the sodium lamp, with the prism turned to the angle of minimum deviation.
- Turn the prism table away by around  $5^\circ$  from the minimum angle so that the ground face of the prism is closer to parallel to the axis of the collimator. In Figure 3, this would mean turning anticlockwise. Now, focus the collimator – adjust the knob on the side of the collimator until the image is a sharp as possible.
- Now turn the prism back through the minimum angle until you are about  $5^\circ$  round on the other side. (Now going clockwise in Figure 3.) Focus the telescope – adjust the knob on the side of the telescope until the image is a sharp as possible.
- Repeat this procedure until no further improvement in the sharpness can be seen.
- The telescope and collimator are now focussed for parallel light.

Original script: Peter Law

Updated script: Peter H Sneddon